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Sapphire Flow

Light Cure Easy Flowable Microhybrid Composite

Sapphire Flow is a visible light cure easy flowing microhybrid composite for anterior and limited posterior restorations and for sealing. It is based on methacrylate-resin and inorganic filler particles of 0.05-1 µm. The total filler load is 61% and the total filler volume 41%. The composite material meets the requirements of DIN EN ISO 4049.

Application of Sapphire Flow

- Class III, IV and V fillings according to Black
- Sealings

Side Effects

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection the cavity floor in cases of deep excavations should be covered with a thin layer of hard setting calciumhydroxid material.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol can impair the hardening of the composite and causes discoloration. Do neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Preparation of Tooth and Color Matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching.

Ascertain the tooth shade while teeth are still moist and select the appropriate Sapphire Flow shade. Sapphire Flow is shaded according to VITA®-Shades.

Application

After isolation, cavity preparation, pulp protection (in case of deep excavations), etching and application of a light cure adhesive system, Sapphire Flow can be applied.

Placement of Sapphire Flow

Apply Sapphire Flow in the selected shade in layers of **no more than 1 mm increments**. Place it by using a non-metallic instrument. Transparent matrix strips may be used.

Curing of Sapphire Flow

Cure for **20 seconds** each layer in case of lighter shades (e.g. А1; А2; С2), cure for **30 seconds** each layer in case of darker shades (e.g. А3.5; В3).

- For minimizing the effects of polymerization shrinkage, curing **in layers** is recommended, if the restoration is **more than 1 mm deep**.

Finishing

Contour with fine diamonds, stones or burs. Polish to high gloss with discs or rubber points. Interproximal finishing is accomplished by fine grit finishing strips.

Storage

Do not store above 25 °C (78 °F). Cool storage prolongs shelf life.

Additional Notes

- Do not use any resin to adjust viscosity of composite restorative material.
- Contact of resin pastes with skin should be avoided, especially by anyone having known resin allergies.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.