

3M™ ESPE™
Stainless Steel Crowns

Unitek™
Stainless Steel Crowns



Crowning Glory

Features and Benefits

- Excellent strength and wear resistance
- Wide array of shapes and sizes
- Comfortable fit
- Made in the U.S.A.
- High quality stainless steel

Consistent, reliable performance for over 40 years.

3M™ ESPE™ Stainless Steel Crowns and Unitek™ Stainless Steel Crowns are #1 in the United States, providing reliable solutions that dentists have trusted for decades.

As a leading innovator in dental products, 3M ESPE offers a full line of prefabricated stainless steel crowns that meet both adult and pediatric patient needs. They are a viable method of single-unit temporization for both short- and long-term coverage.

3M ESPE Stainless Steel Crowns are manufactured with life-like height, contour and occlusal surface. They are pre-primed at the cervical margin.

- Coverage: Short- or long-term
- Dentition: Primary or permanent
- Selection: Molar (primary, permanent)

Unitek Stainless Steel Crowns are prefabricated stainless steel crowns with proven successful clinical use.

- Coverage: Short- or long-term
- Dentition: Primary or permanent
- Selection: Anterior and molar (primary); Bicuspid and molar (permanent)



3M™ ESPE™ Stainless Steel Crowns.

Key Features:

- Pre-trimmed, belled and crimped for fast and easy placement
- Accurately duplicates anatomy for better fit and performance
- Superior longevity compared to multi-surface amalgam fillings
- Permanent molar selection can be used for long-term cases with economic considerations



Unitek™ Stainless Steel Crowns.

Key Features:

- Easy to customize
- Shallow occlusal anatomy means less tooth reduction
- Thick occlusal surface to help prevent bite-through
- Parallel walls allow for operator fit control
- Pre-trimmed to optimum length and contour to save placement time

3M ESPE

Primary					Permanent							
Anterior			Posterior		Posterior							
Centrals Upper and Lower	Laterals Upper and Lower	Cuspids Upper and Lower	Molars		Bicuspid			Molars				
			First Primary Upper and Lower	Second Primary Upper and Lower	First and Second Upper	First Lower	Second Lower	First and Second Upper	First Lower	Second Lower		
					●	●				●	●	
					●	●	●	●	●	●	●	●

Stainless Steel Crown Solutions from 3M ESPE.

3M ESPE and Unitek stainless steel crowns are manufactured in the United States using high quality stainless steel and a proprietary process that results in outstanding physical properties.

The proven performance, fit and esthetics of 3M ESPE and Unitek stainless steel crowns are backed by long-term clinical histories.

For more information on 3M ESPE or Unitek stainless steel crowns, call 1.800.634.2249, visit 3MESPE.com, or contact your 3M ESPE sales representative or 3M ESPE Certified Distributor.

Ordering Information

3M ESPE and Unitek stainless steel crown sets are packaged in a unique, stackable tray that provides convenient, easy storage of prefabricated crowns. Stainless steel crown refills are available and include 5 crowns of one size per box. Crown instruments are also available. Ask your certified 3M ESPE distributor for more information.



3M™ ESPE™ Stainless Steel Primary Molar Crown Set ND96

3M™ ESPE™ Stainless Steel Crowns

Item Number	Product Information
ND96	3M ESPE Stainless Steel Primary Molar Crown Set 96 crowns
P096	3M ESPE Stainless Steel Permanent Molar Crown Set 96 crowns



Unitek™ Stainless Steel Primary Molar Crown Set 902150

Unitek™ Stainless Steel Crowns

Item Number	Product Information
908100	Unitek Stainless Steel Primary Anterior Crown Set 72 crowns
902150	Unitek Stainless Steel Primary Molar Crown Set 112 crowns
902600	Unitek Stainless Steel Permanent Bicuspid Crown Set 84 crowns
902350	Unitek Stainless Steel Permanent Molar Crown Set 84 crowns

Customer Care Center: 1-800-634-2249 www.3MESPE.com



Dental Products

3M Center
Building 275-2SE-03
St. Paul, MN 55144-1000
U.S.A.
1-800-634-2249

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

Please recycle. Printed in U.S.A.
© 3M 2008. All rights reserved.
70-2013-0139-0

3M, ESPE and Unitek are trademarks of 3M or 3M ESPE AG. Used under license in Canada.